



PRINCE MEDICAL

DEVICES



Manufacturer of Medical Devices
**GYNAECOLOGY & MEDICALLY
ASSISTED REPRODUCTION**



PRINCE MEDICAL

DEVICES

3 production sites
in France and Tunisia

10,000 m² of floor space
including 3,000 m² in clean rooms

15 million products manufactured per year



PRINCE MEDICAL, certified NF EN ISO 13485, designs and manufactures single-use medical devices for Gastroenterology, Gynaecology and Medically Assisted Reproduction.

With nearly 30 years of expertise in the healthcare market, the company can rely on technical know-how in the extrusion of thermoplastics and polymers, and has state-of-the-art equipment to meet the needs of markets with demanding applications.

Perfect control of all industrial processes (*moulding, extrusion, production of latex and silicone balloons, assembly of plastic material and steel components, laser cutting, ultrasonic and laser welding, marking, packaging, sterilisation, etc.*) is a guarantee of the high level of quality of PRINCE MEDICAL products.

PRINCE MEDICAL is a subsidiary of the international OMERIN Group and benefits from the dynamism of its 1,700 employees.



The regulatory expertise of PRINCE MEDICAL makes it possible to respond perfectly to current standards and to anticipate their evolution. PRINCE MEDICAL has a long experience in product registration in a number of international markets



PRINCE MEDICAL offers a wide range of sterile Class IIb, Class IIa, Class I sterile and non-sterile medical devices, mainly intended for:

- **Gastroenterology:** balloons catheters, injection needles, cleaning brushes for endoscopes, etc.
- **Gynaecology:** gynaecological forceps (Cheron forceps, Pozzi forceps, etc.), speculums, histological swabbing probes, etc.
- **Medically Assisted Reproduction MAR:** embryo transfer catheters, intrauterine insemination catheter, puncture needle of oocytes, etc.

PRINCE MEDICAL also manufactures medical devices on behalf of its customers for use in urodynamics, anaesthesia, oxygen therapy, as well as external lines for infusion and injection of contrast media in radiology.

The company manufactures and packages its products in a controlled environment and has numerous grey and ISO7 and ISO8 classified clean rooms, guaranteeing a high level of product safety.

The **PRINCE MEDICAL** commercial and technical teams work with passion to put innovation, service and flexibility at the service of its customers.

www.prince-medical.com



OUR PRODUCTS



DIGESTIVE
AND BRONCHIAL
ENDOSCOPY

ENDOLINE[®]
ECHO HOOD[®]
OXYCAP[®]

CLEANING
SWABS AND
BRUSHES

CLEANBRUSH[®]



GYNAECOLOGY
PM-CARE[®]

MEDICALLY ASSISTED
REPRODUCTION
(MAR)

PM-LIFE[®]

GYNAECOLOGY

SAMPLING, SMEAR AND BIOPSY

- PM-CARE® ENDOEXO smear brushPage 8
- PM-CARE® FORBIOPS biopsy forcepsPage 9
- PM-CARE® BRUSH sampling brush.....Page 10
- PM-CARE® ENDOBIOPS histological sampling catheter.....Page 11

INTRA-UTERINE DEVICE (IUD) INSERTION AND REMOVAL

- PM-CARE® EASYCUT string scissors.....Page 13
- PM-CARE® REMOV IUD removal forcepsPage 14
- PM-CARE® INSERT-KIT IUD insertion kitPage 15
- PM-CARE® REMOV-KIT IUD removal kitPage 16
- PM-CARE® RETRIEVAL IUD hook.....Page 17
- PM-CARE® SPECULUM speculum.....Page 24
- PM-CARE® HYSTEROMETRE hystrometer.....Page 25
- PM-CARE® POZZI Pozzi forceps.....Page 26
- PM-CARE® CHERON Cheron forcepsPage 27

HYSTEROGRAPHY AND HYSTEROSONOGRAPHY

- Catheter for hystero- and hysterosonography
1 rigid lumen with foam PM-CARE® 1V-STDPage 19
- Catheter for hystero- and hysterosonography
2 rigid lumens with cup PM-CARE® 2V-CUPPage 20
- Catheter for hystero- and hysterosonography
3 rigid lumens with balloon PM-CARE® 3V-BALLOON.....Page 21
- Hystero- and hysterosonography catheter extension
PM-CARE® EXTENDPage 22

FUNCTIONAL EXPLORATION

- PM-CARE® SPECULUM speculum.....Page 24
- PM-CARE® HYSTEROMETRE hystrometerPage 25
- PM-CARE® CHERON Cheron forcepsPage 26
- PM-CARE® POZZI Pozzi forceps.....Page 27

MAR (MEDICALLY ASSISTED REPRODUCTION)

INTRAUTERINE INSEMINATION

- PM-LIFE® IU embryo transfer setPage 30
- PM-LIFE® PRELEVINJEC aspiration/injection devicePage 31

EMBRYO TRANSFER

- PM-LIFE® FIV IU embryo transfer catheterPage 33
- PM-LIFE® TRANS SET embryo transfer set.....Page 34

OOCYTEPUNCTURE

- PM-LIFE® OVO-KIT oocyte puncture kitPage 36
- PM-LIFE® OVO-PUNC oocyte puncture needle.....Page 37

GYNAECOLOGY

SAMPLING, SMEAR AND BIOPSY	7
INTRA-UTERINE DEVICE (IUD) INSERTION AND REMOVAL	12
HYSTEROGRAPHY AND HYSTEROSONOGRAPHY	18
FUNCTIONAL EXPLORATION	23


SAMPLING, SMEAR AND BIOPSY

- **PM-CARE® ENDOEXO** SMEAR BRUSH 8
- **PM-CARE® FORBIOPS** BIOPSY FORCEPS 9
- **PM-CARE® BRUSH** SAMPLING BRUSH 10
- **PM-CARE® ENDOBIOPS** HISTOLOGICAL SAMPLING CATHETER 11

PM-CARE® ENDOEXO

8

SMEAR BRUSH



The **PM-CARE® ENDOEXO** smear brush is a brush for cervico-vaginal smear enabling simultaneous collection of cells from the endocervix and ectocervix.

ADVANTAGES

- Painless and safe smearing for the patient.
- Long handle for easy smearing and easy manipulation.

AVAILABLE REFERENCE(S)

Reference	Length (cm)
PN3742519	17

ABOUT THIS MEDICAL DEVICE

1 st release date	2013
Single use	Yes
Class	I
Sterile	No
Sterilisation method	None (ethylene oxide decontamination)
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 50 units (two bags of 25)
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com



The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PM-CARE® FORBIOPS**BIOPSY FORCEPS**

9



The PM-CARE® FORBIOPS biopsy forceps is used to perform biopsies on suspect areas of the cervix.

ADVANTAGES

- Scissor-shaped handle for easy manipulation.
- Easy sampling with the jaws of the forceps.
- Detachable pusher for easy retrieval of the sampling.

AVAILABLE REFERENCE(S)

Reference	Length (cm)
PS3600025	26

ABOUT THIS MEDICAL DEVICE

1 st release date	2013
Single use	Yes
Class	Ila
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 20 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PM-CARE® BRUSH

10

SAMPLING BRUSH



**BREAKABLE
TO FACILITATE TRANSPORT
TO THE LABORATORY**

The PM-CARE® BRUSH sampling brush enables cytological swabbing of the endocervix to be performed during a pelvic examination. It is a cytological brush for taking samples, by rotation, for collecting cervical cells in order to diagnose possible pathologies.

ADVANTAGES

- Atraumatic distal end.
- Long handle for easy sampling and easy manipulation.
- Removable part for easy retrieval of the sampling and transport to the laboratory (references PM-CARE® BRUSH breakable).
- Tapered brush configuration for access to narrow openings, for tapered references.

AVAILABLE REFERENCE(S)

Reference	Brush diameter (mm)	Form brush	Sterile	Breakable
PS1643520	3.5	Straight	yes	-
PS3653520	3.5	Straight	yes	yes
PN1643520	3.5	Straight	-	-
PN3653520	3.5	Straight	-	yes
PS1655020	5	Tapered	yes	-
PS3665020	5	Tapered	yes	yes
PN1655020	5	Tapered	-	-
PN3665020	5	Tapered	-	yes

ABOUT THIS MEDICAL DEVICE

1 st release date	2002
Single use	Yes
Class	I
Sterile	*Yes (PS1643520, PS1655020, PS3653520, PS3665020) No (PN1643520, PN1655020, PN3653520, PN3665020)
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 50 units, sterile brushes are packaged in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com



The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PM-CARE® ENDOBIOPS

HISTOLOGICAL SAMPLING CATHETER



The PM-CARE® ENDOBIOPS histological sampling catheter is used to take histological samples from the endometrium.

ADVANTAGES


- Centimetric graduations for better control of catheter insertion.
- Atraumatic tip enabling smooth introduction into the cervical canal.
- Rigid and notched grip tip for a firm grip.
- Efficient aspiration provided by negative pressure.

AVAILABLE REFERENCE(S)

Reference	Model	Diameter (mm / Fr)	Length (cm)
PS1200923	STANDARD	3 / 9	23
PS1210626	MINI	2 / 6	27
PS1220627	CURVED	2 / 6	27
PS3490923	PLUS	3 / 9	23

ABOUT THIS MEDICAL DEVICE


1 st release date	2001
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 units in individual sterile wrapping
Storage	In its original box in a cool and dry place



AVAILABLE MODEL(S)



STANDARD



MINI

- For patients with a narrow cervix or cervical canal
- For post-menopausal patients



CURVED

- For swabbing from the lateral side of the endometrium
- For patients with a non-conventional shaped uterus



PLUS

- To improve swabbing quality
- To take two separate swabbing

GYNAECOLOGY

INTRA-UTERINE DEVICE (IUD) INSERTION AND REMOVAL

- **PM-CARE® EASYCUT** STRING SCISSORS 13
- **PM-CARE® REMOV** IUD REMOVAL FORCEPS 14
- **PM-CARE® INSERT-KIT** IUD INSERTION KIT 15
- **PM-CARE® REMOV-KIT** IUD REMOVAL KIT 16
- **PM-CARE® RETRIEVAL** IUD HOOK 17
- **PM-CARE® SPECULUM** 24
- **PM-CARE® HYSTEROMETRE** 25
- **PM-CARE® POZZI** POZZI FORCEPS 26
- **PM-CARE® CHERON** CHERON FORCEPS 27

PM-CARE® EASYCUT

STRING SCISSORS

13




The PM-CARE® EASYCUT string scissors are used for cutting strings when inserting an intrauterine device. The practitioner can choose to use a pair of lateral or distal string scissors according to their personal preference and ease in the act to be performed.

AVAILABLE REFERENCE(S)

Reference	Model
PS3220022	LATERAL
PS3580022	DISTAL

ABOUT THIS MEDICAL DEVICE

1 st release date	2011
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 20 units in individual sterile wrapping
Storage	In its original box in a cool and dry place



PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.


® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES

PM-CARE® REMOV

14

IUD REMOVAL FORCEPS



The PM-CARE® REMOV IUD removal forceps are used during a gynaecological examination

- When removing an intrauterine device. It holds the strings of the device to pull it down for removal.
- When inserting a new device. It makes it possible to hold a compress for disinfecting and cleaning the work area and/or the cervix soiled by the removal and to prepare the insertion of the new device.

AVAILABLE REFERENCE(S)

Reference	Length (cm)
PS3570025	25

ABOUT THIS MEDICAL DEVICE

1 st release date	2011
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 20 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PM-CARE® INSERT-KIT

IUD INSERTION KIT



CHERON FORCEPS



STRAIGHT POZZI FORCEPS



DISTAL STRING SCISSORS



HYSTEROMETER

The PM-CARE® INSERT-KIT IUD insertion kit is used to insert an intrauterine device in the intrauterine cavity. It is composed of a Pozzi forceps, a Cheron forceps, a pair of string scissors and a hysteroscope.

IUD INSERTION KIT COMPONENTS:

- The straight Pozzi forceps grasps the cervix, performing "traction-fixation" on the cervix to enable the flexion angle of the uterus to be adjusted to 180°. This correction of the flexion angle is necessary in order to insert instruments and fit the intrauterine device without risking uterine perforation.
- The angled Cheron forceps cleans the uterus and the working area.
- Lateral or distal string scissors are used to cut the strings of the intrauterine device. The practitioner can choose to use a pair of lateral or distal string scissors according to their personal preference and ease in the act to be performed.
- The hysteroscope is used to assess the depth of the uterine cavity.

AVAILABLE REFERENCE(S)

Reference	POZZI FORCEPS Model	CHERON FORCEPS Model	STRING SCISSORS Model	HYSTEROMETER Diameter (Fr)
PS4190000	Straight	Angled	Distal	(not supplied)
PS3840000	Straight	Angled	Distal	14
PS3230000	Straight	Angled	Lateral	14

ABOUT THIS MEDICAL DEVICE

1 st release date	2011
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 10 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PM-CARE® REMOV-KIT

IUD REMOVAL KIT



REMOVAL
FORCEPS



STRAIGHT POZZI
FORCEPS



DISTAL STRING
SCISSORS



HYSTEROMETER



The PM-CARE® REMOV-KIT IUD removal kit is used to remove intrauterine devices from the uterine cavity. It is composed of a Pozzi forceps, a removal forceps, a pair of string scissors and a hystrometer.

IUD REMOVAL KIT COMPONENTS:

- The straight Pozzi forceps grasps the cervix, performing "traction-fixation" on the cervix to enable the flexion angle of the uterus to be adjusted to 180°. This correction of the flexion angle is necessary in order to remove the intrauterine device without risking uterine perforation.
- The IUD removal forceps gently remove and cleanse the uterus and working area.
- Distal string scissors are used to cut the strings of the intrauterine device.
- The hystrometer is used to assess the depth of the uterine cavity.

AVAILABLE REFERENCE(S)

Reference	POZZI FORCEPS Model	REMOVAL FORCEPS Furnished	STRING SCISSORS Model	HYSTEROMETER Diameter (Fr)
PS3620000	Straight	Yes	Distal	14

ABOUT THIS MEDICAL DEVICE

1 st release date	2011
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 10 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PM-CARE® RETRIEVAL

IUD HOOK

17



The PM-CARE® RETRIEVAL IUD hook is used to retrieve IUD strings that have been lost in the endocervical canal or the uterine cavity, as well to remove the IUD, if necessary.

AVAILABLE REFERENCE(S)

Reference	Length (mm)
PS3830026	263



ABOUT THIS MEDICAL DEVICE

1 st release date	2014
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 10 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES

GYNAECOLOGY

HYSTEROGRAPHY AND HYSTEROSONOGRAPHY

- **PM-CARE® 1V-STD** HYSTEROGRAPHY AND HYSTEROSONOGRAPHY
1-LUMEN RIGID CATHETER WITH FOAM 19
- **PM-CARE® 2V-STD** HYSTEROGRAPHY AND HYSTEROSONOGRAPHY
2-LUMENS RIGID CATHETER WITH CUP 20
- **PM-CARE® 3V-BALLOON** HYSTEROGRAPHY AND HYSTEROSONOGRAPHY
3-LUMENS MALLEABLE CATHETER WITH BALLOON 21
- **PM-CARE® EXTEND** CATHETER EXTENSION FOR HYSTEROGRAPHY
AND HYSTEROSONOGRAPHY 22

PM-CARE® 1V-STD

1-LUMEN RIGID HYSTEROGRAPHY AND HYSTEROSONOGRAPHY CATHETER WITH FOAM



The PM-CARE® 1V-STD 1-lumen rigid hystero catheter with foam is used to perform hystero-graphy and hystero-sonography and to ensure a tight seal during the examination. The sliding catheter also enables a deeper penetration into the uterine cavity for injection.

ADVANTAGES

- Easy and simple manipulation.
- Sealing ensured by positioning the foam against the cervix.
- Ease of product injection.
- Atraumatic insertion.

AVAILABLE REFERENCE(S)

Reference	Catheter diameter (mm / Fr)	Foam diameter (mm)	Catheter length (mm)	Sliding catheter
PS1870733	2.45 / 7.5	15	330	-
PS2010733	2.45 / 7.5	25	330	-
PS2020733	1.5 / 4.5	15	340	yes

ABOUT THIS MEDICAL DEVICE

1 st release date	2004 2005 for the reference PS2020733
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place



SLIDING

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PM-CARE® 2V-CUP

2-LUMENS RIGID HYSTEROGRAPHY AND HYSTEROSONOGRAPHY CATHETER WITH CUP



The PM-CARE® 2V-CUP 2-lumens rigid hystero catheter with cup is used to perform hystero-graphy and hysterosonography, as well as methylene blue tests.

ADVANTAGES

- The bell-shaped cup is centred on a central core for perfect adhesion of the device, the suction cup function is obtained by applying progressive vacuum.
- The radiologist can use the rigid and malleable manipulator to manipulate the uterus.
- The two lumens can be obturated using a Luer lock valve. The mobile graduated tube enables injections to be performed at the desired location.

AVAILABLE REFERENCE(S)

Reference	Catheter diameter (mm / Fr)	Available cup sizes (mm)	Length (mm)
PS3460845	2.80 / 8.5	23, 28 et 35 (3 sizes provided)	400
PS4292340	2.80 / 8.5	23	400
PS4292840	2.80 / 8.5	28	400
PS4293540	2.80 / 8.5	35	400

ABOUT THIS MEDICAL DEVICE

1 st release date	2012
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PM-CARE® 3V-BALLOON

3-LUMENS MALLEABLE HYSTEROGRAPHY AND HYSTEROSONOGRAPHY CATHETER WITH BALLOON



The PM-CARE® 3V-BALLOON 3-lumens malleable hystero catheter with balloon is used to perform hystero-graphy and hysterosonography. This malleable catheter can be bent and shaped to adapt to the patient's anatomy.

ADVANTAGES

- Easy and simple manipulation.
- Sealing guaranteed by the silicone balloon.
- Standard connector for connecting a syringe for product injection.
- Atraumatic insertion.

AVAILABLE REFERENCE(S)

Reference	Catheter diameter (mm / Fr)	Inflated Balloon diameter (mm)	Balloon insufflation max gas (cc)	Length (mm)
PS3310733	2.66 / 8	13	2.5	230

ABOUT THIS MEDICAL DEVICE

1 st release date	2011
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place



PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PM-CARE® EXTEND

CATHETER EXTENSION FOR HYSTEROGRAPHY AND HYSTEROSONOGRAPHY



The PM-CARE® EXTEND hystero catheter extension connects to catheters for hystero-graphy and hystero-sonography procedures, in order to extend the injection route for easier use by the practitioner.

ADVANTAGES

- Easy and simple manipulation.
- Length adapted for the adjusted use.
- Standard connector for easy and snug connection.

AVAILABLE REFERENCE(S)

Reference	Outer diameter (mm)	Length (cm)
PS0700014	2	140
PS0700002	2	20

ABOUT THIS MEDICAL DEVICE

1 st release date	Ref PS0700014: 2002 Ref PS0700002: 2016
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 10 units in individual sterile wrapping
Storage	In its original box in a cool and dry place



PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com



The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

FUNCTIONAL EXPLORATION

- **PM-CARE® SPECULUM** 24
- **PM-CARE® HYSTEROMETRE** 25
- **PM-CARE® CHERON** CHERON FORCEPS 26
- **PM-CARE® POZZI** POZZI FORCEPS 27

PM-CARE® SPECULUM

24

CUSCO OR COLLIN SPECULUM



The PM-CARE® SPECULUM speculum is used to explore, view, and examine the uterine cervix and vagina, and to access the cervical canal and uterine cavity when performing various examinations and surgical procedures. The Collin speculum enables larger devices to be inserted. The PM-CARE® SPECULUM is made from very high quality medical thermoplastic material.

AVAILABLE REFERENCE(S)

Reference	Model	Size	Diameter (mm)	Colour
PS3592000	Cusco	XS	20	White
PS3592400	Cusco	S	24	Blue
PS3592600	Cusco	M	26	Red
PS3593000	Cusco	L	30	Green
PS3932700	Collin	S	27	White
PS3933000	Collin	M	33	Green

ABOUT THIS MEDICAL DEVICE

1 st release date	Cusco Speculum: 2012 Collin Speculum: 2015
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	• Cusco Speculum: box of 50 units in individual sterile wrapping • Collin Speculum: box of 35 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PM-CARE® HYSTEROMETRE

HYSTEROMETER

25



The PM-CARE® HYSTEROMETRE hystrometer is used to assess uterine depth during gynaecological examinations, in particular for the insertion or removal of a contraceptive device.

AVAILABLE REFERENCE(S)

Reference	Diameter (Fr)	Diameter (mm)	Colour
PS1251000	10	3.3	Blue
PS1251200	12	4	Grey
PS1251400	14	4.6	Green



ABOUT THIS MEDICAL DEVICE

1 st release date	1999
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 or 10 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES

PM-CARE® POZZI

26

POZZI FORCEPS



Pozzi forceps **PM-CARE® POZZI** are used during gynaecological examinations or procedures.

The two main uses:

- Gripping the uterine cervix to hold it still during certain procedures such as cerclage (reduction of the cervical opening), biopsies (removal of cells) or conization (surgical resection of a cone shaped fragment of cervical tissue).
- Provides "traction-fixation" of the uterine cervix in order to adjust the uterine flexion angle to 180°. This correction of the flexion angle is necessary in order to insert instruments and fit or remove the intrauterine device without risking uterine perforation.

PM-CARE® POZZI Pozzi forceps can be straight or angled to adapt to the shape of the uterus. The practitioner can choose to use straight or angled forceps according to their personal preferences and ease in the act to be performed.



AVAILABLE REFERENCE(S)

Reference	Forceps model
PS3180025	Straight
PS3190025	Angled

ABOUT THIS MEDICAL DEVICE

1 st release date	2011
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 20 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PM-CARE® CHERON

CHERON FORCEPS




PM-CARE® CHERON Cheron forceps are used to insert or remove an intrauterine device during gynaecological examinations or procedures. During consultations or gynaecological care, it can be used to grip an aseptic dressing, tampon or compress to apply disinfectant and/or clean the cervix or work area.

PM-CARE® CHERON Cheron forceps can be straight or angled to adapt to the shape of the uterus. The practitioner can choose to use straight or angled forceps according to their personal preferences and ease in the act to be performed.



AVAILABLE REFERENCE(S)

Reference	Forceps model
PS3200025	Straight
PS3210025	Angled



ABOUT THIS MEDICAL DEVICE

1 st release date	2011
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 20 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

MAR

INTRAUTERINE INSEMINATION
EMBRYO TRANSFER
OOCYTE PUNCTURE

29

32

35


INTRAUTERINE INSEMINATION

- **PM-LIFE® IU** INSEMINATION CATHETER 30
- **PM-LIFE® PRELEVINJEC** DEVICE FOR ASPIRATION/INJECTION 31

PM-LIFE® IU

30

INSEMINATION CATHETER



INSERTOR
Stabilises and
guides the catheter

MANDREL
Used to pre-shape the
distal end

The PM-LIFE® IU intrauterine insemination catheter is used to perform insemination during medically assisted reproduction (MAR), by introducing washed and prepared sperm into the uterine cavity. Some product references feature an inserter, to stabilise and guide the catheter during insertion, and/or a mandrel to pre-shape the distal end of the catheter in order to adapt to the patient's morphology.

ADVANTAGES

- Centimetric graduations for better control of catheter insertion.

AVAILABLE REFERENCE(S)

Reference	Diameter for catheter (mm / Fr)	Inner diameter volume (mm / ml)	2 lateral holes diameter (mm)	Useful length (mm)	Length of the distal part (mm)	Inserter	Mandrel
PS1260500	1.65 / 5	1.1 / 0.20	1	175	45	-	-
PS1270500	1.65 / 5	1.1 / 0.20	1	175	45	Yes	-
PS4130500	1.65 / 5	1.1 / 0.20	1	175	45	-	Yes
PS4160500	1.65 / 5	1.1 / 0.20	1	175	45	Yes	Yes
PS1260600	2 / 6	1.6 / 0.40	1.3	175	45	-	-
PS1270600	2 / 6	1.6 / 0.40	1.3	175	45	Yes	-
PS4130600	2 / 6	1.6 / 0.40	1.3	175	45	-	Yes
PS4160600	2 / 6	1.6 / 0.40	1.3	175	45	Yes	Yes
PS2380545	1.6 / 5	1.1 / 0.20	1.1 (distal holes)	175	45	-	Yes

ABOUT THIS MEDICAL DEVICE

1 st release date	2002
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Sterilised via gamma radiation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
DEVICES


The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PM-LIFE® PRELEVINJEC

ASPIRATION/INJECTION DEVICE FOR IUD

31



The PM-LIFE® PRELEVINJEC aspiration/injection device for intrauterine insemination is used to introduce washed and prepared sperm into the uterine cavity without a syringe, as part of artificial intrauterine insemination (IUI) during medically assisted reproduction (MAR).

ADVANTAGES

- Rigid and notched grip tip for a firm grip.
- Atraumatic tip for smooth introduction into the cervical canal.
- Optimal control of injection flow.
- Indicative graduations from 1 to 5 (one graduation corresponding to approximately 100 microlitres).

AVAILABLE REFERENCE(S)

Reference	Inner diameter of the catheter (mm)	Outer diameter of the catheter (mm / Fr)	Inner diameter of the rigid sheath (mm)	Outer diameter of the rigid sheath (mm / Fr)	Length of the catheter (mm)
PS1240623	0.8	2.00 / 6	2.65	3.10 / 9	230



ABOUT THIS MEDICAL DEVICE

1 st release date	2002
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Sterilised via gamma radiation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES

MAR

EMBRYO TRANSFER

- **PM-LIFE® FIV IU** EMBRYO TRANSFER CATHETER 33
- **PM-LIFE® TRANS SET** EMBRYO TRANSFER SET 34

PM-LIFE® FIV IU

EMBRYO TRANSFER CATHETER

33



ECHOGENIC
Echogenic distal part for
better visualisation during
ultrasound



The PM-LIFE® FIV IU embryo transfer catheter is used for intrauterine insemination or to perform a micro-injection (transfer of embryos after in vitro fertilisation). Some product references come with an opaque tube for better ultrasound visualisation, or a mandrel to adapt the catheter to the patient's morphology.

ADVANTAGES

- Centimetric graduations and stopper for better control of catheter insertion.

AVAILABLE REFERENCE(S)

Reference	Catheter diameter (mm / Fr)	Inner diameter / volume (mm / ml)	Excess of the supple catheter (mm)	Useful length (mm)	Tube opaque	Mandrel
PS1280545	1.6 / 5	1.1 / 0.20	45	175	-	-
PS1280555	1.6 / 5	1.1 / 0.20	45	175	Yes	-

ABOUT THIS MEDICAL DEVICE

1 st release date	2002
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Sterilised via gamma radiation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
DEVICES

PM-LIFE® TRANS SET

34

EMBRYO TRANSFER SET



The PM-LIFE® TRANS SET embryo transfer set is used to transfer fertilised embryos "in vitro" into the uterine cavity during a medically assisted reproduction (MAR) procedures.

SET COMPONENTS:

- An insertion catheter. Reference PS3800317 is composed of an "olive" form tip for atraumatic and reliable introduction into the cervix, especially in difficult anatomical conditions.
- An embryo transfer catheter.

AVAILABLE REFERENCE(S)

Reference	INTRODUCTORY CATHETER			TRANSFER CATHETER				
	Diameter (mm / Fr)	Useful length (mm)	Olive	Outer diameter (mm / Fr)	Inner diameter / volume (mm / ml)	Excess of the transfer catheter (mm)	Diameter of the transfer (mm / Fr)	Useful length (mm)
PS1310517	2.0 / 6	160	-	1.5 / 5	1.1 / 0.20	30	1.5 / 5	190
PS1300317	1.6 / 5	175	-	1.0 / 3	0.6 / 0.10	28	1.0 / 3	203
PS3800317	1.6 / 5	175	Yes	1.0 / 3	0.6 / 0.10	57	1.0 / 3	232



OLIVE

Enables atraumatic and reliable insertion

ABOUT THIS MEDICAL DEVICE

1 st release date	2002
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Sterilised via gamma radiation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 10 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.


OOCYTE PUNCTURE

- **PM-LIFE® OVO-KIT** OOCYTE PUNCTURE KIT 36
- **PM-LIFE® OVO-PUNC** OOCYTE PUNCTURE NEEDLE 37

PM-LIFE® OVO-KIT

36

OOCYTE PUNCTURE KIT



The **PM-LIFE® OVO-KIT** oocyte puncture kit is used to perform an ultrasound-guided puncture of the ovarian follicles for the purpose of an in vitro fertilisation. The reference PS4341533 also contains a stopcock (with its male Luer-lock connection base) enabling follicular rinsing.

ADVANTAGES

- The needle is stainless steel with a triple rear bevel and two echogenic areas to capture the ultrasound during ultrasound guidance.

SET COMPONENTS:

- A needle and its attachment base (with tube cover).
- 2 or 3 tubes with their Luer Lock connection bases.
- 1 Falcon-type cap.
- A Stopcock (only for the reference PS4341533).

AVAILABLE REFERENCE(S)

Reference	NEEDLE			
	Inner diameter (mm)	Outer diameter (mm / Gauge)	Useful length (cm)	Stopcock
PS3251533	1.2	1.5 / 17	31.5	-
PS4341533	1.2	1.5 / 17	31.5	Yes

ABOUT THIS MEDICAL DEVICE

1 st release date	2010
Single use	Yes
Class	Ila
Sterile	Yes
Sterilisation method	Sterilised via gamma radiation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	• PS3251533: box of 5 kits in individual sterile wrapping • PS4341533: box of 4 kits in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES



The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PM-LIFE® OVO-PUNC

OOCYTE PUNCTURE NEEDLE

37



The PM-LIFE® OVO-PUNC oocyte puncture needle is used to perform ultrasound-guided puncture of the ovarian follicles for the purpose of an in vitro fertilisation.

ADVANTAGES

- Echogenic distal part for a good visualisation of needle positioning related to the follicle.
- Optimal puncture thanks to the length of the needle and the shape of its bevel.
- A mandrel obturates the needle during insertion to protect oocytes from any contamination

AVAILABLE REFERENCE(S)

Reference	Inner diameter (mm)	Outer diameter (mm / Gauge)	Useful length (cm)
PS1341530	1.2	1.5 / 17	30

ABOUT THIS MEDICAL DEVICE

1 st release date	2002
Single use	Yes
Class	IIa
Sterile	Yes
Sterilisation method	Sterilised via gamma radiation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 10 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE  **MEDICAL**
D E V I C E S

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.





PRINCE  MEDICAL
D E V I C E S

PRINCE MEDICAL SAS
ZA La Sente Du Moulin - 64 rue Gutenberg
60530 Ercuis - France

www.prince-medical.com

Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com