

PRINCE MEDICAL
DEVICES



Manufacturer of Medical Devices

**DIGESTIVE & BRONCHIAL
ENDOSCOPY**

PRINCE MEDICAL

DEVICES

3 production sites
in France and Tunisia

10,000 m² of floor space
including 3,000 m² in clean rooms

15 million products manufactured per year



PRINCE MEDICAL, certified NF EN ISO 13485, designs and manufactures single-use medical devices for Gastroenterology, Gynaecology and Medically Assisted Reproduction.

With nearly 30 years of expertise in the healthcare market, the company can rely on technical know-how in the extrusion of thermoplastics and polymers, and has state-of-the-art equipment to meet the needs of markets with demanding applications.

Perfect control of all industrial processes (*moulding, extrusion, production of latex and silicone balloons, assembly of plastic material and steel components, laser cutting, ultrasonic and laser welding, marking, packaging, sterilisation, etc.*) is a guarantee of the high level of quality of PRINCE MEDICAL products.

PRINCE MEDICAL is a subsidiary of the international OMERIN Group and benefits from the dynamism of its 1,700 employees.



The regulatory expertise of PRINCE MEDICAL makes it possible to respond perfectly to current standards and to anticipate their evolution. PRINCE MEDICAL has a long experience in product registration in a number of international markets



PRINCE MEDICAL offers a wide range of sterile Class IIb, Class IIa, Class I sterile and non-sterile medical devices, mainly intended for:

- **Gastroenterology:** ballons catheters, injection needles, cleaning brushes for endoscopes, etc.
- **Gynaecology:** gynaecological forceps (Cheron forceps, Pozzi forceps, etc.), speculums, histological swabbing probes, etc.
- **Medically Assisted Reproduction MAR:** embryo transfer catheters, intrauterine insemination catheter, puncture needle of oocytes, etc.

PRINCE MEDICAL also manufactures medical devices on behalf of its customers for use in urodynamics, anaesthesia, oxygen therapy, as well as external lines for infusion and injection of contrast media in radiology.

The company manufactures and packages its products in a controlled environment and has numerous grey and ISO7 and ISO8 classified clean rooms, guaranteeing a high level of product safety.

The **PRINCE MEDICAL** commercial and technical teams work with passion to put innovation, service and flexibility at the service of its customers.

www.prince-medical.com





OUR PRODUCTS

DIGESTIVE
AND BRONCHIAL
ENDOSCOPY

ENDOLINE®
ECHO HOOD®
OXYCAP®

CLEANING
SWABS AND
BRUSHES

CLEANBRUSH®



GYNAECOLOGY
PM-CARE®

MEDICALLY
ASSISTED
REPRODUCTION
(MAR)

PM-LIFE®



DIGESTIVE & BRONCHIAL ENDOSCOPY

ERCP Endoscopic retrograde cholangiopancreatography

Balloon catheter

- ENDOLINE® 2V-TRACT 2-lumens balloon catheter Page **12**
- ENDOLINE® 3V-TRACT 3-lumens balloon catheter Page **14**
- ENDOLINE® Y-TRACT Y-clamp fitting..... Page **16**

Nasobiliary kit

- ENDOLINE® NASO-KIT nasobiliary kit..... Page **18**
- ENDOLINE® NASO-DRAIN nasobiliary catheter Page **20**
- ENDOLINE® NASO-CONNECT Luer Lock connector Page **21**
- ENDOLINE® NASO-EXTEND full extension Page **22**
- ENDOLINE® NASO-TUBE nasal tube Page **23**

HAEMOSTASIS

Injection needle

- ENDOLINE® 1V-JECT 1-way injection needle..... Page **26**
- ENDOLINE® 2V-JECT 2-ways injection needle..... Page **28**

Polypectomy

- ENDOLINE® LOOP polypectomy snare Page **30**
- ENDOLINE® COLD-LOOP cold snare..... Page **32**
- ENDOLINE® SPRAY catheter spray..... Page **33**

CYTOLOGICAL SAMPLING

- ENDOLINE® CYTO-BRUSH cytology brush Page **36**

GASTRIC BALLOON EXTRACTION

- ENDOLINE® PUNC drainage needle..... Page **40**
- ENDOLINE® VIPER extraction forceps Page **41**

ECHO-ENDOSCOPY

- ECHO HOOD® APERTURE-1 balloon with 1 opening Page **44**
- ECHO HOOD® APERTURE-2 balloon with 2 openings Page **45**

FOREIGN BODIES

- ENDOLINE® HOOD hood..... Page **48**

ACCESSORIES

- ENDOLINE® VALVE valve for operator channel Page **50**
- OXYCAP® EVO mouthpiece Page **52**

CLEANING SWABS AND BRUSHES

- CLEANBRUSH® ENDO cleaning brush for endoscope operator channelPage 54
- CLEANBRUSH® ENDO-XL cleaning brush for rigid instruments Page 56
- CLEANBRUSH® PISTON cleaning brush for endoscope valvesPage 57
- CLEANBRUSH® EXO-B brush for external cleaning of endoscopesPage 58
- CLEANBRUSH® EXO-S sponge for external cleaning of endoscopes.....Page 59
- CLEANBRUSH® KIT-2 cleaning swab kit - 2 partsPage 60
- CLEANBRUSH® KIT-3 cleaning swab kit - 3 partsPage 62
- CLEANBRUSH® KIT-4 cleaning swab kit - 4 parts.....Page 64
- CLEANBRUSH® KIT-5 cleaning swab kit - 5 partsPage 64
- CLEANBRUSH® TRACH brush for tracheal cannulaPage 66





DIGESTIVE & BRONCHIAL ENDOSCOPY

ERCP	11
HAEMOSTASIS	25
CYTOLOGICAL SAMPLING	35
GASTRIC BALLOON EXTRACTION	39
ECHO-ENDOSCOPY	43
FOREIGN BODIES	47
ACCESSORIES	49



DIGESTIVE & BRONCHIAL ENDOSCOPY

ERCP

• ENDOLINE® 2V-TRACT 2-LUMENS BALLOON CATHETER	12
• ENDOLINE® 3V-TRACT 3-LUMENS BALLOON CATHETER	14
• ENDOLINE® Y-TRACT Y-CLAMP FITTING	16
• ENDOLINE® NASO-KIT NASOBILIARY KIT	18
• ENDOLINE® NASO-DRAIN NASOBILIARY CATHETER	20
• ENDOLINE® NASO-CONNECT LUER-LOCK CONNECTOR	21
• ENDOLINE® NASO-EXTEND FULL EXTENSION	22
• ENDOLINE® NASO-TUBE NASAL TUBE	23

ENDOLINE® 2V-TRACT

12

2-LUMENS BALLOON CATHETER



QUICK EXCHANGE SYSTEM
Makes it easier to insert catheter into echo-endoscope
AS2340720

WELDED DERIVATIONS

BONDED DERIVATIONS

TRANSPARENT TUBE
AS4460720

The ENDOLINE® 2V-TRACT 2-lumens balloon catheter enables the removal of gallstones and injection of contrast material or any other medicated substance into the bile ducts. This medical device is intended mainly for persons with hepatic colic following the appearance of gallstones.

ADVANTAGES

- Product compatible for use with the main endoscopes on the market.
- Radiopaque.
- Standard connector for connecting a syringe for contrast medium injection.
- Presence of a stainless steel ring and stainless steel ligatures on the balloon reinforcing the radiopacity of the distal part for better balloon localisation.

ADDITIONAL ACCESSORIES

The balloon catheter can be used with the ENDOLINE® Y-TRACT Y-clamp fitting. Provides practitioners with 2 injection ports on balloon catheters with a single injection port. It also enables the fixation of the guide wire after catheter placement.

AVAILABLE BALLOONS



LATEX BALLOON CATHETER

PRECAUTIONS FOR USE
DO NOT USE IN PATIENTS ALLERGIC TO LATEX



SILICONE BALLOON CATHETER

for patients allergic to latex.



ABOUT THIS MEDICAL DEVICE

1 st release date	1996
Single use	Yes
Class	Ila
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Latex present	Yes except reference AS0920720
Lifetime	2 years except reference AS0920720 5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com



The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® 2V-TRACT

LATEX
BALLOON

Reference	Catheter diameter (mm/Fr)	Inflated balloon diameter (mm)	Balloon inflation max gas (mL)	Balloon injection max liquid (mL)	Wire-guide recommended (inch)	Syringe (cc)	Useful length (cm)	Option
> Bonded derivations								
AS0030520	1.70 / 5	12	3.0	1.5	.025	3	200	-
AS1100720	2.35 / 7 1.70 / 5	16	4.5	2.5	.035	5	200	-
AS1910720	2.35 / 7	12	3.0	1.5	.035	3	200	-
AS0030720	2.35 / 7	16	4.5	2.5	.035	5	200	-
AS4460720	2.35 / 7	16	4.5	2.5	.035	5	200	Transparent
AS0280720	2.35 / 7	18	5.5	3.5	.035	5	200	-
> Welded derivations								
AS2340720	2.35 / 7	16	4.5	2.5	.035	5	200	Quick exchange system
AS2040720	2.35 / 7 1.70 / 5	16	4.5	2.5	.035	5	200	Tube cover
AS2130720	2.35 / 7	16	4.5	2.5	.035	5	200	Tube cover

SILICONE
BALLOON

> Bonded derivations								
AS0920720	2.35 / 7	16	4.5	2.5	.035	5	200	-

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® 3V-TRACT

14

3-LUMENS BALLOON CATHETER



TRANSPARENT
TUBE
AS4560720

The ENDOLINE® 3V-TRACT 3-lumens balloon catheter enables the removal of gallstones and injection of contrast material or any other medicated substance into the bile ducts. This medical device is intended mainly for persons with hepatic colic following the appearance of gallstones.

ADVANTAGES

- Product compatible for use with the main endoscopes on the market.
- Radiopaque.
- Standard connector for connecting a syringe for contrast medium injection.
- Stainless steel ring and stainless steel ligatures on the balloon enhance radiopacity of the distal part for better balloon identification.
- The guide path and the contrast path are different.

AVAILABLE BALLOONS



LATEX BALLOON CATHETER

PRECAUTIONS FOR USE
DO NOT USE IN PATIENTS
ALLERGIC TO LATEX



SILICONE BALLOON CATHETER

for patients allergic
to latex.



ABOUT THIS MEDICAL DEVICE

1 st release date	1996
Single use	Yes
Class	Ila
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Latex present	Yes except reference AS3860720
Lifetime	2 years except reference AS3860720 5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® 3V-TRACT

LATEX
BALLOON

Reference	Catheter diameter (mm/Fr)	Inflated balloon diameter (mm)	Balloon inflation max gas (mL)	Balloon injection max liquid (mL)	Wire-guide recommended (inch)	Syringe (cc)	Useful length (cm)	Option
AS1420720	2.35 / 7	16	4.5	2.5	.035	5	200	-
AS4320720	2.35 / 7	16	4.5	2.5	.035	5	200	Transparent
AS4320725	2.35 / 7	16	4.5	2.5	.035	5	200	Transparent
AS4560720	2.35 / 7 1.70 / 5	16	4.5	2.5	.035	5	200	Transparent

SILICONE
BALLOON

AS3860720	2.35 / 7	16	4.5	2.5	.035	5	200	-
-----------	----------	----	-----	-----	------	---	-----	---

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com


The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® Y-TRACT

16

Y CLAMP FITTING



The ENDOLINE® Y-TRACT Y-clamp fitting provides practitioners with 2 injection ports on a balloon catheter with 1 injection port, and also enables the fixation of the guide wire after catheter placement.

AVAILABLE REFERENCE(S)

Reference	Description
AS1130300	CH3 Y-clamp fitting

ABOUT THIS MEDICAL DEVICE

1 st release date	1998
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 50 units in individual sterile wrapping
Storage	In its original box in a cool and dry place



PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.



® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.



ENDOLINE® NASO-KIT

18

NASOBILIARY KIT



The ENDOLINE® NASO-KIT nasobiliary kit enables the aspiration of bile or pancreatic juices. This medical device is intended primarily for persons with obstructive jaundice, acute cholangitis or extrahepatic cholestasis.

NASOBILIARY KIT COMPONENTS (components can be sold separately)

- A nasobiliary catheter for aspirating bile juice.
- A Luer-Lock connector mounted on the nasobiliary catheter enabling biliary or pancreatic juice aspiration. It connects the nasobiliary catheter to the full extension.
- A full extension to drain bile or pancreatic secretions in an external receiving bag or to rinse the catheter as needed.
- A nasal tube to retrieve the nasobiliary catheter in the larynx and to remove it through the nose.

ADVANTAGES

- Axial and lateral perforations in the distal part allow for efficient drainage of the hepatobiliary tract.
- Easy to position catheter in the bile ducts thanks to the radiopaque distal end (with gold ring).

AVAILABLE FORMATS (NASOBILIARY CATHETER)



CURVE 1



CURVE 2



CURVE 6

ABOUT THIS MEDICAL DEVICE

1 st release date	1996
Single use	Yes
Class	Ila
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 1 unit in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® NASO-KIT

Reference	NASOBILIARY CATHETER		CONNECTOR LUER-LOCK	NASAL TUBE		EXTENSION
	Catheter diameter (Fr)	Catheter length (cm)	Connector Luer-Lock	Nasal tube Int. diam. / Ext. diam. (mm / mm)	Length of the nasal tube (cm)	Length of the extension (cm)
CURVE 1						
AS0540500	5	250	grey	2.0 / 3.0	50	135
AS0540700	7	250	orange	2.6 / 3.7	50	135
AS0541000	10	250	white	3.7 / 4.9	50	135
CURVE 2						
AS0560500	5	250	grey	2,0 / 3,0	50	135
AS0560700	7	250	orange	2,6 / 3,7	50	135
CURVE 6						
AS0580500	5	250	grey	2,0 / 3,0	50	135
AS0580700	7	250	orange	2,6 / 3,7	50	135

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com


The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® NASO-DRAIN

20

NASOBILIARY CATHETER



The ENDOLINE® NASO-DRAIN nasobiliary catheter is used to aspirate bile or pancreatic juices. This medical device is intended primarily for persons with obstructive jaundice, acute cholangitis or extrahepatic cholestasis.

ADVANTAGES

- Axial and lateral perforations in the distal part enable efficient drainage of the hepatobiliary tract.
- Easy positioning of the catheter in the bile ducts thanks to the radiopaque distal extremity (with gold ring).

AVAILABLE FORMS



CURVE 1



CURVE 2



CURVE 6

AVAILABLE REFERENCE(S)

	Reference	Catheter diameter (Fr)	Useful length (cm)
CURVE 1	AS0530525	5	250
	AS0530725	7	250
	AS0531025	10	250
CURVE 2	AS0550525	5	250
	AS0550725	7	250
CURVE 6	AS0570525	5	250
	AS0570725	7	250

ABOUT THIS MEDICAL DEVICE

1 st release date	1996
Single use	Yes
Class	Ila
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES


The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® NASO-CONNECT

LUER-LOCK CONNECTOR

21



The ENDOLINE® NASO-CONNECT Luer-Lock connector enables the aspiration of bile or pancreatic juices. It also allows to remove the endoscope while leaving the catheter in place. This component is part of the ENDOLINE® NASO-KIT nasobiliary kit (for more information on the kit, please consult the data sheet page 18).

AVAILABLE REFERENCE(S)

Reference	Diameter (Fr)	Colour
AS0050500	5	grey
AS0050600	6	blue
AS0050700	7	orange
AS0051000	10	white



ABOUT THIS MEDICAL DEVICE

1 st release date	1997
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 50 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.


® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES

ENDOLINE® NASO-EXTEND

22

FULL EXTENSION




The ENDOLINE® NASO-EXTEND full extension drains bile or pancreatic secretions into an external receiving bag, and can also be used to rinse the nasobiliary catheter. This single use medical device is part of the ENDOLINE® NASO-KIT nasobiliary kit (for more information on the kit, please consult the data sheet page 18).

It has a 3-ways stopcock:

- 1 way to the patient, which is connected to the proximal end of the nasobiliary catheter via the Luer-Lock connector.
- 1 way for rinsing of the nasobiliary catheter, when required.
- 1 way to collect bile in an external receiving bag.

AVAILABLE REFERENCE(S)

Reference	Length (cm)
AS0210013	135



ABOUT THIS MEDICAL DEVICE

1 st release date	1997
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES


The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® NASO-TUBE

NASAL TUBE


23



The ENDOLINE® NASO-TUBE nasal tube is used to retrieve a nasobiliary catheter from the larynx and remove it through the nose. This component is part of the ENDOLINE® NASO-KIT nasobiliary kit (for more information on the kit, please consult the data sheet page 18).

AVAILABLE REFERENCE(S)

Reference	Nasal tube Int. diam. / Ext. diam. (mm / mm)	Length of the nasal tube (cm)	Compatible with nasobiliary catheter diameter (Fr)
ASO220550	2.0 / 3.0	50	5
ASO220750	2.6 / 3.7	50	7
ASO221050	3.7 / 4.9	50	10



ABOUT THIS MEDICAL DEVICE

1 st release date	1996
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 10 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES



DIGESTIVE & BRONCHIAL ENDOSCOPY

25


HAEMOSTASIS

- **ENDOLINE® 1V-JECT** 1-WAY INJECTION NEEDLE 26
- **ENDOLINE® 2V-JECT** 2-WAYS INJECTION NEEDLE 28
- **ENDOLINE® LOOP** POLYPECTOMY SNARE 30
- **ENDOLINE® COLD-LOOP** COLD SNARE 32
- **ENDOLINE® SPRAY** CATHETER SPRAY 33

ENDOLINE® 1V - JECT

26

1-WAY INJECTION NEEDLE



PLASTIC TIP



METAL TIP



The ENDOLINE® 1V-JECT 1-way injection needle has been designed specifically and exclusively for haemostasis of gastrointestinal haemorrhages, sclerosis of oesophageal varices and inflating and marking polyps.

The needle is only suitable for use with flexible endoscopes.


ADVANTAGES

- Radiopaque distal part visible on the screen, for models with a metal tip.
- Ergonomically shaped handle for easy manipulation.
- Product compatible for use with the main endoscopes on the market.
- Luer-Lock connector for connecting a syringe for product injection


AVAILABLE HANDLES



1.0
The most widely used handle in the world, used and approved for over 20 years!



2.0
Evolution of the handle 1.0, thicker for a better grip in the hand.



3.0
Optimal ergonomics and lighter weight, this handle has been designed with the user and the environment in mind!

ABOUT THIS MEDICAL DEVICE

1 st release date	1996
Single use	Yes
Class	Ila
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Erceuil - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com



The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® 1V - JECT

HANDLE 1.0

> Metal Tip

Reference	Catheter diameter (mm/Fr)	Needle diameter (mm)	(G)	Useful length (cm)	Needle outlet length (mm)
AS1990718	2,35/7	0,7	22	180	5
AS1990723	2,35/7	0,7	22	230	5
AS1991018	2,35/7	1	19	180	5,5
AS1991023	2,35/7	1	19	230	5,5

HANDLE 2.0

> Plastic Tip

AS2940518	1,75/5,3	0,5	25	180	4
AS2940523	1,75/5,3	0,5	25	230	4
AS2940718	2,35/7	0,7	22	180	5
AS2940723	2,35/7	0,7	22	230	5
AS2941018	2,35/7	1	19	180	5,5
AS2941023	2,35/7	1	19	230	5,5

> Metal Tip

AS3010718	2,35/7	0,7	22	180	5
AS3010723	2,35/7	0,7	22	230	5
AS3011018	2,35/7	1	19	180	5,5
AS3011023	2,35/7	1	19	230	5,5
AS3720523	2,35/7	0,5	25	230	3
AS4310523	2,35/7	0,5	25	230	4,5

HANDLE 3.0


> Plastic Tip

AS4480518	1,75/5,3	0,5	25	180	4
AS4480523	1,75/5,3	0,5	25	230	4
AS4480718	2,35/7	0,7	22	180	5
AS4480723	2,35/7	0,7	22	230	5
AS4481018	2,35/7	1	19	180	5
AS4481023	2,35/7	1	19	230	5

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

ENDOLINE® 2V - JECT


2 WAYS INJECTION NEEDLE



PLASTIC TIP



METAL TIP



The ENDOLINE® 2V-JECT 2-ways injection needle has been designed specifically and exclusively for haemostasis of gastrointestinal haemorrhages, sclerosis of oesophageal varices and inflating and marking polyps.

The needle is only suitable for use with flexible endoscopes.

The 2nd way can be used to clean cleaning and/or rinse the procedure area.

ADVANTAGES

- Radiopaque distal part visible on the screen, for models with a metal tip.
- Ergonomically shaped handle for easy manipulation.
- Product compatible for use with the main endoscopes on the market.
- Luer-Lock connector for connecting a syringe for product injection.

AVAILABLE HANDLES



1.0
The most widely used handle in the world, used and approved for over 20 years!



2.0
Evolution of the handle 1.0, thicker for a better grip in the hand.

ABOUT THIS MEDICAL DEVICE

1 st release date	1996
Single use	Yes
Class	IIa
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Erceuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

ENDOLINE® 2V - JECT

HANDLE 1.0

> Metal Tip

Reference	Catheter diameter (mm/Fr)	Needle diameter (mm) (G)		Useful length (cm)	Length of the needle outlet (mm)
AS1840718	2,35 / 7	0,7	22	180	5
AS1840723	2,35 / 7	0,7	22	230	5
AS1841018	2,35 / 7	1	19	180	5,5
AS1841023	2,35 / 7	1	19	230	5,5

HANDLE 2.0

> Plastic Tip

AS2830518	1,75/5,3	0,5	25	180	4
AS2830523	1,75/5,3	0,5	25	230	4
AS2830718	2,35/7	0,7	22	180	5
AS2830723	2,35/7	0,7	22	230	5
AS2830728	2,35/7	0,7	22	280	5
AS2831018	2,35/7	1	19	180	5,5
AS2831023	2,35/7	1	19	230	5,5

> Metal Tip


AS2900718	2,35/7	0,7	22	180	5
AS2900723	2,35/7	0,7	22	230	5
AS2900728	2,35/7	0,7	22	280	5
AS2901018	2,35/7	1	19	180	5,5
AS2901023	2,35/7	1	19	230	5,5

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

ENDOLINE® LOOP

30

POLYPECTOMY SNARE



**INNER CABLE
AVAILABLE IN 0.67 MM AND
0.9 MM DIAMETER**

**The use of a larger diameter
inner cable
(0.9 mm) enables the distal part
to be heated to a higher temper-
ature*

The ENDOLINE® LOOP polypectomy snare enables the resection of sessile/flat or pedicled polyps (≥ 6 mm) in the digestive tract via an upper or lower endoscopy. This snare is connected to the diathermy energy generator via the insert plug on the handle's movable slider.

ADVANTAGES

- Three fingers grip for easy manipulation.
- Radiopaque distal part for perfect visibility.
- Product compatible for use with the main endoscopes on the market.
- Connector for connection to a diathermic energy generator.

AVAILABLE FORMS



LARGE OVAL
Snare with flexible or rigid cable
Ø 30 x 50 mm



SMALL OVAL
Snare with flexible or rigid cable
Ø 15 x 25 mm



MONOWIRE OVAL
Snare with flexible
cable
Ø 25 x 50 mm
or 17 x 45 mm



HEXAGONAL
Snare with flexible or rigid cable
Ø 30 x 50 mm



ASYMMETRIC
Snare with flexible cable
Width, 15 mm or 25 mm

ABOUT THIS MEDICAL DEVICE

1 st release date	2001
Single use	Yes
Class	IIB
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® LOOP

LARGE OVAL

Reference	Type of snare	Diameter of the inner cable (mm)	Shape diameter (cable/rod) (mm)	Cable / Rod	Sheath diameter (mm)	Sheath length (cm)	Opening (mm)
AS2865023	Flexible	0,67	0,39	Cable	2,35	230	30 x 50
AS3755023	Rigid	0,9	0,45	Cable	2,35	230	30 x 50

SMALL OVAL

AS2922523	Flexible	0,67	0,39	Cable	2,35	230	15 x 25
AS3762523	Rigid	0,9	0,45	Cable	2,35	230	15 x 25

MONOWIRE OVAL

AS1961723	Flexible	0,67	0,25	Rod	2,35	230	17 x 45
AS1972523	Flexible	0,67	0,30	Rod	2,35	230	25 x 50

HEXAGONAL

AS2985023	Flexible	0,67	0,39	Cable	2,35	230	30 x 50
AS3775023	Rigid	0,9	0,45	Cable	2,35	230	30 x 50

ASYMMETRIC

AS3265023	Flexible	0,67	0,39	Cable	2,35	230	25
AS3961523	Flexible	0,67	0,39	Cable	2,35	230	15

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

ENDOLINE® COLD-LOOP

32

COLD SNARE



The ENDOLINE® COLD-LOOP cold snare is intended for endoscopic resection of small sessile/flat (≤ 9 mm) polyps present in the gastrointestinal tract via an upper or lower endoscopy without use of diathermic energy.

ADVANTAGES

- Three fingers grip for easy manipulation.
- Radiopaque distal part for perfect visibility.
- Product compatible for use with the main endoscopes on the market.

AVAILABLE REFERENCE(S)

Reference	Sheath diameter (mm)	Sheath length (cm)	Opening (mm)
AS4372523	2.35	230	10 x 25

ABOUT THIS MEDICAL DEVICE

1 st release date	2016
Single use	Yes
Class	IIa
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® SPRAY

CATHETER SPRAY

33



The ENDOLINE® SPRAY catheter spray enables staining products to be injected into the digestive system's mucous membrane. This injection is used for mucosal analysis in case of lesions or polyps in the digestive system, in patients with a relevant medical history or for diagnosis purposes.

ADVANTAGES


- Easy one-handed operation.
- Product compatible for use with the main endoscopes on the market.
- Luer-Lock connector for connecting a syringe for injection of staining products.

AVAILABLE REFERENCE(S)

Reference	Catheter diameter (mm)	Useful length (cm)
AS2791823	1.8	230

ABOUT THIS MEDICAL DEVICE

1 st release date	2009
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place



PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES



DIGESTIVE & BRONCHIAL ENDOSCOPY

35

CYTOLOGICAL SAMPLING


- **ENDOLINE® CYTO-BRUSH** CYTOLOGY BRUSH

36

ENDOLINE® CYTO-BRUSH

36

CYTOLOGY BRUSH



The ENDOLINE® CYTO-BRUSH cytology brush enables cytological samples to be taken from the mucous membranes of the digestive tract, bronchial tubes and bile ducts for the diagnosis of pathologies such as infectious diseases or cancer.

ADVANTAGES

- Handle with thumb ring for easy manipulation.
- Radiopaque distal part for perfect visibility.
- Includes scissors (depending on product reference).
- Product compatible for use with the main endoscopes on the market.

ADDITIONAL ACCESSORIES

Scissors to cut the brush at the end of the examination and place it in a storage bottle for transport to the test lab.

The references concerned: AS3673022, AS3673016, AS3672512, AS3672522



ABOUT THIS MEDICAL DEVICE

1 st release date	1996
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® CYTO-BRUSH

WITHOUT ACCESSORIES

Reference (cm)	Sheath length (mm)	Sheath diameter (mm)	Brush diameter
AS0373022	220	2.35	3
AS0373016	160	2.35	3
AS0372512	120	1.75	2.5
AS0372522	220	1.75	2.5

WITH ACCESSORIES
(SCISSORS)

AS3673022	220	2.35	3
AS3673016	160	2.35	3
AS3672512	120	1.75	2.5
AS3672522	220	1.75	2.5

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.



DIGESTIVE & BRONCHIAL ENDOSCOPY

39


GASTRIC BALLOON EXTRACTION

- **ENDOLINE® PUNC** DRAINAGE NEEDLE 40
- **ENDOLINE® VIPER** EXTRACTION FORCEPS 41

ENDOLINE® PUNC

40

DRAINAGE NEEDLE



The ENDOLINE® PUNC drainage needle is used to drain gastric balloons.

ADVANTAGES

- Ergonomically shaped handle for easy manipulation.
- Radiopaque distal part for perfect visibility.
- Luer-Lock connector for connecting a syringe to drain air and water from the gastric balloons.
- Product compatible for use with the main endoscopes on the market.

AVAILABLE REFERENCE(S)

Reference	Catheter diameter (mm / Fr)	Tip diameter (mm)	Useful length (cm)
AS2640718	2.33 / 7	1	180

AVAILABLE AS A KIT WITH EXTRACTION FORCEPS

ENDOLINE® PUNC&VIPER / *Extraction kit*

Reference: AS2710718

ABOUT THIS MEDICAL DEVICE

1 st release date	2008
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES


The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® VIPER

EXTRACTION FORCEPS

41



The ENDOLINE® VIPER double-grasp extraction forceps is used for the extraction of emptied intragastric balloons.

ADVANTAGES

- Three fingers grip for easy manipulation.
- Radiopaque distal part for perfect visibility.
- Product compatible for use with the main endoscopes on the market.
- Reinforced with a sheath spring to support balloon resistance during extraction.

AVAILABLE REFERENCE(S)


Reference	Sheath diameter (mm / Fr)	Sheath length (cm)
AS2290718	2.3/7	180

AVAILABLE AS A KIT WITH THE DRAINAGE NEEDLE**ENDOLINE® PUNC&VIPER / Extraction kit**

Reference: AS2710718

ABOUT THIS MEDICAL DEVICE

1 st release date	2008
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place



PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES



DIGESTIVE & BRONCHIAL ENDOSCOPY

43

ECHO-ENDOSCOPY

- **ECHO HOOD® APERTURE-1** BALLOON 1 OPENING 44
- **ECHO HOOD® APERTURE-2** BALLOON 2 OPENINGS 45

ECHO HOOD® APERTURE-1

44

BALLOON 1 OPENING



The ECHO HOOD® APERTURE-1 balloon for echo-endoscope fits on the ultrasound probe to make contact with the gastrointestinal wall and thus provide a good ultrasound image.

ADVANTAGES

- Can withstand handling.
- Adequate thickness for the passage of ultrasound waves.
- Product compatible for use with the main endoscopes on the market.

AVAILABLE REFERENCE(S)

Reference	Compatible endoscopes	Balloon length (mm)	Opening size (mm)
AS1470000	ECHO-DUODENOSCOPE Pentax	34	8.7
AS1830000	ECHO-DUODENOSCOPE Olympus	15	6.4
AS3680000	ECHO-DUODENOSCOPE Olympus 180	15	6.5
AS2800000	ECHO-ENDOSCOPE Fujinon convex	18	5.8
AS3560000	ECHO-ENDOSCOPE Fujinon convex	18	6



ABOUT THIS MEDICAL DEVICE

1 st release date	1998
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	Yes
Lifetime	2 years
Certification	CE
Packaging	Box of 20 units in individual sterile wrapping
Storage	In its original box, protected from ultraviolet rays, in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES


The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ECHO HOOD® APERTURE-2

BALLOON 2 OPENINGS

45



The ECHO HOOD® APERTURE-2 balloon for echo-endoscope fits on the ultrasound probe to make contact with the gastrointestinal wall and thus provide a good ultrasound image.

ADVANTAGES

- Can withstand handling.
- Adequate thickness for the passage of ultrasound waves.
- Product compatible for use with the main endoscopes on the market.

PRECAUTIONS FOR USE

DO NOT USE IN PATIENTS ALLERGIC TO LATEX

AVAILABLE REFERENCE(S)

Reference	Compatible endoscopes	Length of the balloon (mm)	Size of 1 st opening (mm)	Size of 2 nd opening (mm)
AS1070000	ECHO-DUODENOSCOPE Olympus	16	5.4	2.7
AS2070000	ECHO-DUODENOSCOPE Pentax	19	6.5	5
AS3710000	ECHO-DUODENOSCOPE Pentax	22	1.8	4.2
AS2660000	ECHO-ENDOSCOPE Fujinon radial	19	5.8	5.8
AS3610000	ECHO-ENDOSCOPE Fujinon radial	19	6	6
AS2230000	ECHO-ENDOSCOPE series 160 Olympus	17	5.4	2.7
AS1460000	ECHOCOLOSCOPE Olympus	17	6	8.8
AS3280000	ECHO-BRONCHOSCOPE Olympus	12	1.3	3.3
AS3940000	ECHO-BRONCHOSCOPE Pentax	13.5	3.1	1.15

ABOUT THIS MEDICAL DEVICE

1 st release date	1998
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	Yes
Lifetime	2 years
Certification	CE
Packaging	Box of 20 units in individual sterile wrapping
Storage	In its original box, protected from sunlight, in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES



DIGESTIVE & BRONCHIAL ENDOSCOPY

47

FOREIGN BODIES

- **ENDOLINE® HOOD HOOD**

48

ENDOLINE® HOOD

48

HOOD



The ENDOLINE® HOOD fits on the distal end of most endoscopes and protects the cardia, oesophagus and pharynx during endoscopic foreign body removal. It also prevents the foreign body from escaping.

ADVANTAGES

- Bell shape to completely cover the foreign body.
- Easy fixation and compatible with main endoscopes on the market.
- Secured by the attached elastic band.

PRECAUTIONS FOR USE

DO NOT USE IN PATIENTS ALLERGIC TO LATEX

AVAILABLE REFERENCE(S)

Reference	Material	Distal tip diameter (mm)	Proximal tip diameter (mm)	Useful length (mm)
AS0420000	Latex	40	8.3	75

ABOUT THIS MEDICAL DEVICE

1 st release date	1998
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	Yes
Lifetime	2 years
Certification	CE
Packaging	Box of 10 units in individual sterile wrapping
Storage	In its original box, protected from ultraviolet rays, in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

DIGESTIVE & BRONCHIAL ENDOSCOPY


ACCESSORIES

- **ENDOLINE® VALVE** VALVE FOR WORKING CHANNEL 50
- **OXYCAP® EVO** MOUTHPIECE 52

ENDOLINE® VALVE

50

WORKING CHANNEL VALVE



VALVE WITH LATERAL IRRIGATION FITTING



VALVE WITH LATERAL IRRIGATION FITTING AND NON-RETURN VALVE

The ENDOLINE® VALVE working channel valve is a phosphorescent valve that fits on the endoscope's working channel. It gives accessories access to the working channel, guarantees a tight seal during the examination and eliminates the risk of cross-contamination between patients.

ADVANTAGES

- Narrow closure to prevent leakage.
- Product compatible for use with the main endoscopes on the market.
- References with trans-valve fitting connecting to commonly used syringe models.
- References with lateral irrigation fitting for lavage through the working channel while maintaining the accessory in place.
- References with integrated non-return valve to prevent gastric reflux in the endoscope.

ABOUT THIS MEDICAL DEVICE

1 st release date	2007
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 50 or 100 units in individual sterile wrapping, depending on product reference
Storage	In its original box in a cool and dry place



VALVE ONLY



VALVE WITH FITTING

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com



The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ECHO HOOD® VALVE

DIGESTIVE & BRONCHIAL

Reference	Compatible endoscopes	Supplied with a connector	Colour	Packaging
AS2430000	Olympus, Fujinon	-	orange	100
AS2440000	Pentax	-	yellow	100
AS2500000	Olympus, Fujinon	trans-valve	orange	100
AS2510000	Pentax	trans-valve	yellow	100
AS2450000	Olympus, Fujinon	lateral irrigation	orange	50
AS2460000	Pentax	lateral irrigation	yellow	50
AS3810000	Olympus, Fujinon	lateral irrigation & AR valve	orange	50
AS3820000	Pentax	lateral irrigation AR & valve	yellow	50

BRONCHIAL ENDOSCOPY

AS2690000	Olympus, Fujinon	-	blue	100
AS3020000	Pentax	-	pink	100
AS3350000	Olympus, Fujinon	trans-valve	blue	100
AS3400000	Pentax	trans-valve	pink	100

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com


The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

OXYCAP® EVO

52

MOUTHPIECE



MOUTHPIECE
WITH 2 CONNECTORS
MOUNTED



MOUTHPIECE
WITH BAND
AND 2 CONNECTORS
MOUNTED

The OXYCAP® EVO mouthpiece is used in upper digestive endoscopies (gastroscopy, duodenoscopy, etc.) to protect the patient's teeth, as well as to protect endoscopes passing through the oral cavity. The connectors allow continuous oxygenation during the procedure.

MOUTHPIECES

Reference	Configuration	Packaging
AS4200000	Only	75
AS4210000	With band	75
AS4250000	With 2 connectors mounted	75
AS4260000	With band and 2 connectors mounted	75

ADDITIONAL ACCESSORIES

Reference	Configuration	Packaging
AS4220000	Band only	75
AS4270000	Connector with Luer	10
AS4280000	Connector with filter	10
AS4230000	White connector only	100
AS4390000	Yellow connector only	100

ABOUT THIS MEDICAL DEVICE

1 st release date	2017
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	2 years
Certification	CE
Packaging	Box of 10, 75 or 100 units in individual sterile wrapping depending on product reference
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

DIGESTIVE & BRONCHIAL ENDOSCOPY

CLEANING SWABS AND BRUSHES

- **CLEANBRUSH® ENDO**
CLEANING BRUSH FOR ENDOSCOPE WORKING CHANNEL 54
- **CLEANBRUSH® ENDO-XL**
CLEANING BRUSH FOR RIGID INSTRUMENTS 56
- **CLEANBRUSH® PISTON**
CLEANING BRUSH FOR ENDOSCOPE VALVES 57
- **CLEANBRUSH® EXO-B**
BRUSH FOR EXTERNAL ENDOSCOPE CLEANING 58
- **CLEANBRUSH® EXO-S**
SPONGE FOR EXTERNAL ENDOSCOPE CLEANING 59
- **CLEANBRUSH® KIT-2**
2-PART CLEANING SWAB KIT 60
- **CLEANBRUSH® KIT-3**
3-PART CLEANING SWAB KIT 62
- **CLEANBRUSH® KIT-4**
4-PART CLEANING SWAB KIT 64
- **CLEANBRUSH® KIT-5**
5-ELEMENT CLEANING SWAB KIT 64
- **CLEANBRUSH® TRACH**
BRUSH FOR TRACHEAL CANNULA 66

CLEANBRUSH® ENDO

54




CLEANING BRUSH FOR ENDOSCOPE WORKING CHANNEL



CLEANBRUSH® ENDO cleaning brushes are used to clean the

- > working channels, air/water channels in flexible endoscopes,
- > working channels, air/water channels in rigid endoscopes,
- > O-type duodenoscope elevator.

They are suitable for use with the main gastroscopes, colonoscopes, duodenoscopes, enteroscopes, bronchoscopes, cystoscopes, ureteroscopes, nephroscopes, hysteroscopes and naso-laryngoscopes on the market.



ABOUT THIS MEDICAL DEVICE

1 st release date	1999
Single use	Yes
Class	I
Sterile	No (ethylene oxide decontamination)
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 50 units in non-sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

AVAILABLE REFERENCE(S)

> CLEANING THE WORKING CHANNEL ON FLEXIBLE ENDOSCOPES

Reference	Brush diameter (mm)	Number of brushes	Sheath diameter (mm)	Recommended diameter for working channel of the endoscope	Length (cm)
AN3272011	2	1	0.8	From 1 mm to 1.6 mm	110
AN3272023	2	1	0.8	From 1 mm to 1.6 mm	230
AN1050106	1.8	1	1	From 1.2 mm to 1.6 mm	60
AN1050110	1.8	1	1	From 1.2 mm to 1.6 mm	100
AN1050123	1.8	1	1	From 1.2 mm to 1.6 mm	230
AN1050135	1.8	1	1	From 1.2 mm to 1.6 mm	350
AN1050138	1.8	1	1	From 1.2 mm to 1.6 mm	380
AN1060110	1.8	2	1	From 1.2 mm to 1.6 mm	100
AN1060123	1.8	2	1	From 1.2 mm to 1.6 mm	230
AN1053006	3	1	1.60	From 1.8 mm to 2.6 mm	60
AN1053007	3	1	1.60	From 1.8 mm to 2.6 mm	75
AN1053010	3	1	1.60	From 1.8 mm to 2.6 mm	100
AN1053018	3	1	1.60	From 1.8 mm to 2.6 mm	180
AN1053023	3	1	1.60	From 1.8 mm to 2.6 mm	230
AN1053028	3	1	1.60	From 1.8 mm to 2.6 mm	280
AN1063010	3	2	1.60	From 1.8 mm to 2.6 mm	100
AN1063018	3	2	1.60	From 1.8 mm to 2.6 mm	180
AN1063023	3	2	1.60	From 1.8 mm to 2.6 mm	230
AN1063028	3	2	1.60	From 1.8 mm to 2.6 mm	280
AN2083018	11& 3	2	1.60	From 1.8 mm to 2.6 mm	180
AN2083023	11& 3	2	1.60	From 1.8 mm to 2.6 mm	230
AN1055002	5	1	1.80	From 2.7 mm to 4.6 mm	20
AN1055006	5	1	1.80	From 2.7 mm to 4.6 mm	60

> CLEANING THE WORKING CHANNEL ON RIGID ENDOSCOPES

Reference	Brush Diameter (mm)	Number of brushes	Sheath diameter (mm)	Recommended diameter for working channel of the endoscope	Length (cm)
AN1053004	3	1	1.60	From 1.8 mm to 2.6 mm	40
AN1055004	5	1	1.80	From 2.7 mm to 4.6 mm	40
AN1051104	11	1	1.80	From 6.1 mm to 10.6 mm	40

Reference	Brush diameter (mm)	Number of brushes	Sheath diameter (mm)	Recommended diameter for working channel of the endoscope	Length (cm)
AN1055010	5	1	1.80	From 2.7 mm to 4.6 mm	100
AN1055018	5	1	1.80	From 2.7 mm to 4.6 mm	180
AN1055020	5	1	1.80	From 2.7 mm to 4.6 mm	200
AN1055023	5	1	1.80	From 2.7 mm to 4.6 mm	230
AN1055028	5	1	1.80	From 2.7 mm to 4.6 mm	280
AN1055030	5	1	1.80	From 2.7 mm to 4.6 mm	300
AN1065010	5	2	1.80	From 2.7 mm to 4.6 mm	100
AN1065018	5	2	1.80	From 2.7 mm to 4.6 mm	180
AN1065020	5	2	1.80	From 2.7 mm to 4.6 mm	200
AN1065023	5	2	1.80	From 2.7 mm to 4.6 mm	230
AN1065028	5	2	1.80	From 2.7 mm to 4.6 mm	280
AN1065030	5	2	1.80	From 2.7 mm to 4.6 mm	300
AN2085023	11 & 5	2	1.80	From 2.7 mm to 4.6 mm	230
AN1785023	5	3	1.80	From 2.7 mm to 4.6 mm	230
AN3525023	5	2	1.80	From 2.7 mm to 4.6 mm	230
AN1056506	6.5	1	1.80	From 4.7 mm to 6 mm	60
AN1056523	6.5	1	1.80	From 4.7 mm to 6 mm	230
AN1066523	6.5	2	1.80	From 4.7 mm to 6 mm	230
AN1068023	8	2	1.80	From 4.2 mm to 7.6 mm	230
AN1068028	8	2	1.80	From 4.2 mm to 7.6 mm	280

> OTHER CLEANING

Reference	Brush diameter (mm)	Number of brushes	Type of cleaning	Length (cm)
AN4302000	2	1	0-type duodenoscope elevator	6



PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

CLEANBRUSH® ENDO-XL

56

CLEANING BRUSH FOR RIGID INSTRUMENTS



The **CLEANBRUSH® ENDO-XL** cleaning brushes are intended for cleaning the channels of rigid instruments before sterilisation, such as surgical trocars.

They are also suitable for cleaning rigid endoscopes used in urology, gynaecology and ENT.

AVAILABLE REFERENCE(S)

Reference	Brush diameter (mm)	Number of brushes	Sheath diameter (mm)	Recommended diameter for channel to be cleaned	Length (cm)	Packaging
AN3435506	5.5	1	3.3	From 4.1 to 5 mm	60	20
AN3438006	8	1	3.3	From 5.1 mm to 7 mm	60	20
AN3431106	11	1	3.3	From 7.1 to 9 mm	60	20
AN3431406	14	1	3.3	From 9.1 mm to 12 mm	60	10
AN3431806	18	1	3.3	From 12.1 to 16 mm	60	10
AN2575006	5 11	2	3.3	From 2.7 to 4.5 mm From 6 to 10 mm	60	10



ABOUT THIS MEDICAL DEVICE

1 st release date	1999
Single use	Yes
Class	I
Sterile	No (ethylene oxide decontamination)
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 10 or 20 units in non-sterile wrapping according to the reference
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES


The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

CLEANBRUSH® PISTON

CLEANING BRUSH FOR ENDOSCOPE VALVES

57



The CLEANBRUSH® PISTON cleaning brushes are intended for cleaning the valve and valve cages of flexible endoscopes.

They are suitable for use with the main gastroscopes, colonoscopes, duodenoscopes, enteroscopes, bronchoscopes, cystoscopes, ureteroscopes, nephroscopes, hysteroscopes and naso-laryngoscopes on the market.

AVAILABLE REFERENCE(S)

Reference	Brush diameter (mm)	Number of brushes	Type of cleaning	Length (cm)
AN1061110	11 & 5	2	For valves and valve cages	15.5



ABOUT THIS MEDICAL DEVICE

1 st release date	1999
Single use	Yes
Class	I
Sterile	No (ethylene oxide decontamination)
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 50 units in non-sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES

CLEANBRUSH® EXO-B

58

BRUSH FOR EXTERNAL ENDOSCOPE CLEANING



The CLEANBRUSH® EXO-B external brushes are intended for cleaning the external surface of endoscopes.

AVAILABLE REFERENCE(S)

Reference	Bristles length (cm)	Useful Length (cm)	Brush shape	Bristle type
AN2580000	1	2 x 8	U-shaped	Supple
AN3952200	2	5.5	Straight with handle	Supple
AN3953500	1.5	5.5	Straight with handle	Rigid

ABOUT THIS MEDICAL DEVICE

1 st release date	1999
Single use	Yes
Class	I
Sterile	No (ethylene oxide decontamination)
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 units in non-sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com


PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

CLEANBRUSH® EXO-S

SPONGE FOR EXTERNAL ENDOSCOPE CLEANING



The CLEANBRUSH® EXO-S sponge is intended for cleaning the external surface of endoscopes.

AVAILABLE REFERENCE(S)

Reference	Available width (cm)	Available length (cm)
AN1790000	5	13



ABOUT THIS MEDICAL DEVICE

1 st release date	1999
Single use	Yes
Class	I
Sterile	No (ethylene oxide decontamination)
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 units in non-sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.


® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.



CLEANBRUSH® KIT-2

60

2-PART CLEANING SWAB KIT



The CLEANBRUSH® KIT-2 swab kits are intended for cleaning the working channels, valves and valve cages of flexible endoscopes.

They are suitable for use with most major brands of gastroscopes, colonoscopes, duodenoscopes, enteroscopes, bronchoscopes, cystoscopes, ureteroscopes, nephroscopes, hysteroscopes and naso-laryngoscopes.

KIT COMPONENTS (SEE THE TABLE FOR THE DETAILED COMPONENTS BY REFERENCE)

- CLEANBRUSH® ENDO: cleaning brush for working channel
- CLEANBRUSH® PISTON: cleaning brush for valves, valve cages
- CLEANBRUSH® EXO-B: brush for external cleaning
- CLEANBRUSH® EXO- S: sponge for external cleaning

ABOUT THIS MEDICAL DEVICE

1 st release date	1999
Single use	Yes
Class	I
Sterile	No (ethylene oxide decontamination)
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 or 50 units depending on product reference
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

CLEANBRUSH® KIT-2

61

REFERENCES FOR CLEANING BRUSHES

AVAILABLE REFERENCE(S)

Kit reference	For working channel	For valves	Packaging
AN1850110	AN1050110 (with 1 brush Ø 1.8 mm L100 cm)		
AN1853010	AN1053010 (with 1 brush Ø 3 mm L100 cm)		
AN1853018	AN1053018 (with 1 brush Ø 3 mm L180 cm)		
AN1855010	AN1055010 (with 1 brush Ø 5 mm L100 cm)		
AN1855018	AN1055018 (with 1 brush Ø 5 mm L180 cm)		
AN1855020	AN1055020 (with 1 brush Ø 5 mm L200 cm)		
AN1855023	AN1055023 (with 1 brush Ø 5 mm L230 cm)		
AN1856523	AN1056523 (with 1 brush Ø 6.5 mm L230 cm)		
AN1860110	AN1060110 (with 2 brushes Ø 1.8 mm L100 cm)		
AN1863010	AN1063010 (with 2 brushes Ø 3 mm L100 cm)		
AN1863018	AN1063018 (with 2 brushes Ø 3 mm L180 cm)		
AN1865010	AN1065010 (with 2 brushes Ø 5 mm L100 cm)		
AN1865018	AN1065018 (with 2 brushes Ø 5 mm L180 cm)		
AN1865020	AN1065020 (with 2 brushes Ø 5 mm L200 cm)		
AN1865023	AN1065023 (with 2 brushes Ø 5 mm L230 cm)		
AN1866523	AN1066523 (with 2 brushes Ø 6.5 mm L230 cm)		
		AN1061110 (with 2 brushes Ø11 & 5mm, L15.5 cm)	50

KIT COMPONENT REFERENCES

Kit reference	Brush for external cleaning	Sponge for external cleaning	Packaging
AN2750000	AN2580000 (U-shaped with soft hair, useful length 2 x 8 cm)	AN1790000	25

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.


® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES

CLEANBRUSH® KIT-3

62

3-PART CLEANING SWAB KIT



The CLEANBRUSH® KIT-3 swab kits are intended for cleaning the working channels, valves and valve cages of flexible endoscopes.

They are suitable for use with the main gastroscopes, colonoscopes, duodenoscopes, enteroscopes, bronchoscopes, cystoscopes, ureteroscopes, nephroscopes, hysteroscopes and naso-laryngoscopes on the market.

KIT COMPONENTS (SEE THE TABLE FOR THE DETAILED COMPONENTS BY REFERENCE)

- CLEANBRUSH® ENDO: cleaning brush for working channel
- CLEANBRUSH® ENDO (AN1055002, AN1056523): cleaning brush for aspiration channel
- CLEANBRUSH® ENDO (AN1050123): cleaning brush for air/water channel, and/or Water Jet and/or balloon inflation channel
- CLEANBRUSH® PISTON: cleaning brush for valves, valve cages
- CLEANBRUSH® EXO-B: brush for external cleaning

ABOUT THIS MEDICAL DEVICE

1 st release date	1999
Single use	Yes
Class	I
Sterile	No (ethylene oxide decontamination)
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 or 50 units depending on product reference
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

CLEANBRUSH® KIT-3

AVAILABLE REFERENCE(S)

REFERENCES FOR CLEANING BRUSHES

OTHER KIT COMPONENTS

Kit reference	For working channel	For aspiration channel	For valves	For air/water channel, water jet, balloon inflation channel	External cleaning brush	Packaging	
AN1883018	AN1053018 (with 1 brush Ø 3 mm L180 cm)					50	
AN1885023	AN1055023 (with 1 brush Ø 5 mm L230 cm)						
AN1886523	AN1056523 (with 1 brush Ø 6.5 mm L230 cm)						
AN1893018	AN1063018 (with 2 brushes Ø 3 mm L180 cm)						
AN1895023	AN1065023 (with 2 brushes Ø 5 mm L230 cm)						
AN1896523	AN1066523 (with 2 brushes Ø 6.5 mm L230 cm)						
AN2770110	AN1050110 (with 1 brush Ø 1.8 mm L100 cm)						
AN2780110	AN1060110 (with 2 brushes Ø 1.8 mm L100 cm)						
AN2183010	AN1053010 (with 1 brush Ø 3 mm L100 cm)						AN1055002 (with 1 brush Ø 5 mm, L20 cm)
AN2183018	AN1053018 (with 1 brush Ø 3 mm L180 cm)						
AN2193010	AN1063010 (with 2 brushes Ø 3 mm L100 cm)				50		
AN2193018	AN1063018 (with 2 brushes Ø 3 mm L180 cm)						
AN2953018	AN1053018 (with 1 brush Ø 3 mm L180 cm)						
AN2963018	AN1063018 (with 2 brushes Ø 3 mm L180 cm)						
AN4475018	AN1055018 (with 1 brush Ø 5 mm L180 cm)						
AN4475023	AN1055023 (with 1 brush Ø 5 mm L230 cm)						
AN4495018	AN1065018 (with 2 brushes Ø 5 mm L180 cm)						
AN4495023	AN1065023 (with 2 brushes Ø 5 mm L230 cm)						
AN2520110	AN1050110 (with 1 brush Ø 1.8 mm L100 cm)						
AN2523007	AN1053007 (with 1 brush Ø 3 mm L75 cm)						
AN2523010	AN1053010 (with 1 brush Ø 3 mm L100 cm)				25		
AN2523018	AN1053018 (with 1 brush Ø 3 mm L180 cm)						
AN2525010	AN1055010 (with 1 brush Ø 5 mm L100 cm)						
AN2525023	AN1055023 (with 1 brush Ø 5 mm L230 cm)						
AN2525030	AN1055030 (with 1 brush Ø 5 mm L300 cm)						
AN2530110	AN1060110 (with 2 brushes Ø 1.8 mm L100 cm)						
AN2533010	AN1063010 (with 2 brushes Ø 3 mm L100 cm)						
AN2533018	AN1063018 (with 2 brushes Ø 3 mm L180 cm)						
AN2535010	AN1065010 (with 2 brushes Ø 5 mm L100 cm)						
AN2535023	AN1065023 (with 2 brushes Ø 5 mm L230 cm)						
AN2535030	AN1065030 (with 2 brushes Ø 5 mm L300 cm)						


The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

CLEANBRUSH® KIT-4, KIT-5

64

CLEANING SWAB KIT – 4 OR 5 PARTS



The CLEANBRUSH® KIT-4, KIT-5 swab kits are intended for cleaning the working channels, valves and valve cages of flexible endoscopes.

They are suitable for use with the main gastroscopes, colonoscopes, duodenoscopes, enteroscopes, bronchoscopes, cystoscopes, ureteroscopes, nephroscopes, hysteroscopes and naso-laryngoscopes on the market.

KIT COMPONENTS (SEE THE TABLE FOR THE DETAILED COMPONENTS BY REFERENCE)

- CLEANBRUSH® ENDO: cleaning brush for working channel
- CLEANBRUSH® ENDO (AN1055002, AN1055023, AN1056523): cleaning brush for the aspiration channel
- CLEANBRUSH® ENDO (AN1050123, AN1050138): cleaning brush for air/water channel, and/or Water Jet and/or balloon inflation channel
- CLEANBRUSH® PISTON: cleaning brush for valves, valve cages
- CLEANBRUSH® EXO-B: brush for external cleaning

ABOUT THIS MEDICAL DEVICE

1 st release date	1999
Single use	Yes
Class	I
Sterile	No (ethylene oxide decontamination)
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 or 50 units depending on product reference
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

CLEANBRUSH® KIT-4, KIT-5

AVAILABLE REFERENCE(S)

> KIT with 4 parts

Kit reference	REFERENCES FOR CLEANING BRUSHES				OTHER KIT COMPONENTS	Packaging
	For working channel	For aspiration channel	For valves	For air/water channel, water jet, balloon inflation channel	External cleaning brush	
AN1903018	AN1053018 (with 1 brush Ø 3 mm L180 cm)	AN1055023 (with 1 brush Ø 5 mm L230 cm)				50
AN1983018	AN1063018 (with 2 brushes Ø 3 mm L180 cm)			AN1050123 (with 1 brush Ø 1.8 mm L230 cm)		
AN4535018	AN1065018 (with 2 brushes Ø 5 mm L180 cm)	AN1056523 (with 1 brush Ø 6.5 mm L230 cm)				25
AN2543010	AN1053010 (with 1 brush Ø 3 mm L100 cm)	AN1055002 (with 1 brush Ø 5 mm L20 cm)	AN1061110 (with 2 brushes Ø 11 & 5 mm, L15.5 cm)			
AN2553010	AN1063010 (with 2 brushes Ø 3 mm L100 cm)					
AN4505023	AN1055023 (with 1 brush Ø 5 mm L230 cm)					
AN4515023	AN1065023 (with 2 brushes Ø 5 mm L230 cm)	AN1056523 (with 1 brush Ø 6.5 mm L230 cm)			AN2580000 (U-shaped with soft hair, useful length 2 x 8 cm)	
AN4505018	AN1055018 (with 1 brush Ø 5 mm L180 cm)					
AN4515018	AN1065018 (with 2 brushes Ø 5 mm L180 cm)					
AN2543018	AN1053018 (with 1 brush Ø 3 mm L180 cm)					
AN2545023	AN1055023 (with 1 brush Ø 5 mm L230 cm)			AN1050123 (with 1 brush Ø 1.8 mm L230 cm)		
AN3335023	AN1065023 (with 2 brushes Ø 5 mm L230 cm)					
AN2995023	AN1055023 (with 1 brush Ø 5 mm L230 cm)			AN1050138 (with 1 brush Ø 1.8 mm L380 cm)		

AVAILABLE REFERENCE(S)

> KIT with 5 parts

Kit reference	REFERENCES FOR CLEANING BRUSHES				OTHER KIT COMPONENTS	Packaging
	For working channel	For aspiration channel	For valves	For air/water channel, water jet, balloon inflation channel	External cleaning brush	
AN4545018	AN1055018 (with 1 brush Ø 5 mm L180 cm)					25
AN4555018	AN1065018 (with 2 brushes Ø 5 mm L180 cm)	AN1056523 (with 1 brush Ø 6.5 mm L230 cm)	AN1061110 (with 2 brushes Ø 11 & 5 mm, L15.5 cm)	AN1050123 (with 1 brush Ø 1.8 mm L230 cm)	AN2580000 (U-shaped with soft hair, useful length 2 x 8 cm)	

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

CLEANBRUSH® TRACH

66


BRUSH FOR TRACHEAL CANNULA



The CLEANBRUSH® TRACH tracheal cannula brushes are intended for cleaning tracheostomy cannulas.

AVAILABLE REFERENCE(S)

Reference	Brush diameter (mm)	Sheath diameter (mm)	Length (cm)	Packaging
AS2201002	10	3.3	23	25
AS2201402	14	3.3	23	25
AS2201802	18	3.3	23	20



ABOUT THIS MEDICAL DEVICE

1 st release date	2007
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 20 or 25 units in sterile wrapping depending on product reference
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE  MEDICAL
D E V I C E S

PRINCE MEDICAL SAS
ZA La Sente Du Moulin - 64 rue Gutenberg
60530 Ercuis - France

www.prince-medical.com

Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com